

instructions: MITO PRO+ FLOOR STAND

www.mitoredlight.com

© 2020 Mito Red Light, Inc. All rights reserved.



ABOUT US

Located in Scottsdale, Arizona, **Mito Red Light**[®] is dedicated to bringing the highest quality red light therapy lights to market at the best value.

Having extensively studied and realized the benefits of red light therapy firsthand, Mito Red Light was brought into existence in order to provide consumers truly affordable, yet highly innovative red light therapy options.

If you have any further questions, please feel free to contact us at: **info@mitoredlight.com**

You can also call our phone number Monday - Friday 9am - 5pm MST +1866-861-6486

Mail correspondence can be sent to: Mito Red Light, Inc 9319 N 94th Way, Suite 400 Scottsdale, AZ 85258

Instructions: MitoPRO+ Floor Stand

Step 1: Unbox your stand and unpackage the provided hardware.



Step 2: Thread the provided casters into the base of the floor stand and hand tighten. (optional).



Step 3: Flip the stand onto its casters and install the 2 side braces and tighten using the provided screws with a Philips head screwdriver. (Ensure the holes line up as shown in the picture).





Instructions: MitoPRO+ Floor Stand - Single Panel

Step 1: You will need to remove the lower two rubber feet from the back of your unit using a small Philips head screwdriver before installation into the floor stand. **Step 2:** Insert your light into the Floor stand and ensure it is seated properly.



Important note: ensure the panel is fully seated into the bottom of the stand by pushing downwards on the top of the panel.





‼ mitoredlight

www.mitoredlight.com



User Manual and Instructional Videos This user manual and instructional videos can be found here: https://mitoredlight.com/pages/user-manual

Follow Us!









light /mitore

/mitoredlightofficial

@MitoRedLight

@mitoredlightofficial